

# LAPAROSCOPY INSTRUMENTS






**LAPARO MED**









PASSION FOR HEALING

[www.laparomed.com](http://www.laparomed.com)




# INTRODUCING ELBOWLESS HAND INSTRUMENTS





## ATRAUMATIC GRASPER





## DISSECTING GRASPING FORCEP



## SCISSOR




## TRAUMATIC GRASPER




## VERESS NEEDLE



Veress Needle  
Dia 2.1mm, Length 12cm  
■ VN-21-120



Dia 2.1mm, Length 14cm  
■ VN-21-140



Veress Needle  
Dia 2.9mm, Length 14cm  
■ VN-29-140

## TROCAR CANNULA

LAPAROMED TROCAR OFFERS THE ULTIMATE FACILITY, FLEXIBILITY, VALUE AND QUICK CUSTOMIZATION FOR ANY SURGICAL PROCEDURE



Trocar Cannula CFT  
Spiral Pyramidal  
Dia 6mm, Length 9cm  
■ CFT-06-90-10-PC



Trocar Cannula CFT  
Plain Pyramidal  
Dia 10 mm, Length 10cm  
■ CFT-10-100-20-PC



Trocar Cannula  
Spiral Plain  
Dia 5mm, Length 9cm  
■ SFT-05-90-20-PC



Trocar Cannula CFT  
Plain Pyramidal  
Dia 7mm, Length 9cm  
■ SFT-07-90-20-PC



Trocar Cannula  
Plain  
Dia 10mm, Length 10cm  
■ SFT-10-100-20-PC

## SAFETY TROCAR



Safety Trocar VN Type CFT  
Dia 5mm, Length 16.5cm  
■ STC-06-165-VN



Safety Trocar VN Type CFT  
Dia 10mm, Length 20cm  
■ STC-10-200-VN




Safety Trocar with Lock CFT  
Dia 10mm, Length 20cm  
■ STC-10-200




Hassan Rod & Cone  
For CFT 10cm  
■ THC-10-100-BC

## REDUCER




Reducer Metal Sleeve  
10/5mm (for CFT)  
■ RCM-105-175



Reducer Metal Sleeve  
10/7mm (for CFT)  
■ RCM-107-175



Reducer Metal Sleeve  
7/5mm (for CFT)  
■ RCM-075-165



Reducer Metal Sleeve  
05/3mm (for CFT)  
■ RCM-053-155



Reducer Flip On  
10/5mm  
■ RUF-105-20



Reducer Flip On  
10/7mm  
■ RUF-107-20



Reducer Flip On  
07/5mm  
■ RUF-075-20




Reducer Flip On  
05/3mm  
■ RUF-053-20




# DISSECTING /GRASPING FORCEPS



Maryland Small Dissecting & Grasping Forcep FHP  
Dia 5mm, Length 30cm  
■ **DMS-05-30-01**




Maryland Medium Dissecting & Grasping Forcep FHP  
Dia 5mm, Length 30cm  
■ **DMM-05-30-01**




Maryland Medium Diamond Dissecting & Grasping Forcep FHP  
Dia 5mm, Length 30cm  
■ **DMMD-05-30-01**




Maryland Long Dissecting & Grasping Forcep FHP  
Dia 5mm, Length 30cm  
■ **DML-05-30-01**



Maryland Medium Dissecting & Grasping Forcep FHP  
Dia 10mm, Length 30cm  
■ **DMM-10-30-01**




Maryland 45 ° Dissecting & Grasping Forcep FHP  
Dia 5mm, Length 30cm  
■ **DMA-05-30-45°-01**




Mixer 90° Dissecting & Grasping Forcep FHP  
Dia 5mm, Length 30cm  
■ **DMR-05-30-90°-01**




Mixer 90° Dissecting & Grasping Forcep FHP  
Dia 10mm, Length 30cm  
■ **DMR-10-30-90°-01**



Dolphin Dissecting Grasping & Forcep FHP  
Dia 5mm, Length 30cm  
■ **DDO-05-30-01**



Duckbill Dissecting & Grasping Forcep FHP  
Dia 5mm, Length 30cm  
■ **DDB-05-30-01**



full range of control tip dissectors for enhance precision used in laparoscopic surgery.


# SCISSORS

sharp scissors for surgery to help optimize cutting performance and reusable scissors.




Micro Straight Scissor (SA) FHP  
Dia 5mm, Length 30cm

■ SMS-05-30-01-SA



Straight Scissor (SA) FHP  
Dia 5mm, Length 30cm

■ SSS-05-30-01-SA




Straight Scissor (DA) FHP  
Dia 5mm, Length 30cm

■ SSD-05-30-01



Curve Scissor Metzanbaum  
FHP Dia 5mm, 30 Length

■ SCM-05-30-01



Curve Scissor Maryland Type  
FHP Dia 5mm, Length 30cm

■ SCD-05-30-01



Curve Scissor Short FHP  
Dia 5mm, Length 30cm

■ SCS-05-30-01




Hook Scissor (SA) FHP  
Dia 5mm, Length 30cm

■ SHS-05-30-01-SA




Cardia Achalasia Scissor (SA) FHP  
Dia 10mm, Length 30cm


■ SCA-10-30-01-SA




sharp scissors for surgery to help optimize cutting performance and reusable scissors.




Micro Straight Scissor (SA) FHP  
Dia 5mm, Length 30cm  
■ SMS-05-30-01-SA




Straight Scissor (SA) FHP  
Dia 5mm, Length 30cm  
■ SSS-05-30-01-SA



Straight Scissor (DA) FHP  
Dia 5mm, Length 30cm  
■ SSD-05-30-01



Curve Scissor Metzbaum  
FHP Dia 5mm, 30 Length  
■ SCM-05-30-01



Curve Scissor Maryland Type  
FHP Dia 5mm, Length 30cm  
■ SCD-05-30-01




Curve Scissor Short FHP  
Dia 5mm, Length 30cm  
■ SCS-05-30-01



Hook Scissor (SA) FHP  
Dia 5mm, Length 30cm  
■ SHS-05-30-01-SA




Cardia Achalasia Scissor (SA) FHP  
Dia 10mm, Length 30cm  
■ SCA-10-30-01-SA



# ATRAUMATIC GRASPERS




Plain Long Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **APL-05-30-03**




Allies Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AAS-05-30-03**




Alligator Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AGA-05-30-03**



Endoclinch Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AGE-05-30-03**




Universal Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AGU-05-30-03**




Fenestrated Small (SA)  
Grasping Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AFS-05-30-03**



Fenestrated Medium (DA)  
Grasping Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AFM-05-30-03**




Fenestrated Long Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AFL-05-30-03**




Babcock Atrogrip Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **ABA-05-30-03**



Babcock Plate Grasping  
Forcep FHBR  
Dia 10mm, Length 30cm  
■ **ABP-10-30-03**



Oval Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AOS-05-30-03**



Tube Holder Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **ATH-05-30-03**




Endowave Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AGW-05-30-03**



Fenestrated Double  
Slit Straight (SA) Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AFDS-05-30-03-SA**



Fenestrated Double  
Slit Curve (SA) Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **AFD-05-30-03-SA**



Babcock Paddle Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ **ABP-05-30-03**



## TRAUMATIC GRASPERS



2x4 Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ T24-05-30-03



2x3 Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ T23-05-30-03



Davinchi Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ TGD-05-30-03




Tenaculam (SA) Grasping  
Forcep FHBR  
Dia 5mm, Length 30cm  
■ T11-05-30-03-SA


## BIOPSY FORCEPS




Biopsy Grasping Forcep  
Spoon FHP  
Dia 5mm, Length 30cm  
■ BFS-05-30-01



Biopsy Grasping Forcep  
(Spoon with Spike) FHP  
Dia 5mm, Length 30cm  
■ BSS-05-30-01



Clamp Intestinal  
Straight FHBR  
Dia 5mm, Length 30cm  
■ CIS-05-30-03



Clamp Intestinal  
Curved FHBR  
Dia 5mm, Length 30cm  
■ CIC-05-30-03


## INTESTINAL CLAMP

Detachable range of efficient ergonomically designed Traumatic graspers to

## GALL BLADDER EXTRACTOR




Claw Forceps Traumatic Dia 10mm, Length 30cm with Metal Handle Spring Loaded  
■ GBE-10-30-M1C



## STONE REMOVAL



Stone Removal Forceps Dia 10mm, Length 30cm, with Metal Handle Spring Loaded  
■ SRF-10-30-M1C



# CLIP APPLICATOR




Clip Applicator Dia 5mm,  
Length 30cm, Medium LT200

**C2-05-30-08-E**



Clip Applicator 10mm, 30cm  
Medium Large LT300


**C3-10-30-07-E**



L M/L

Clip Applicator 2-in-1  
Dia 10mm, Length 30cm  
for LT 300/400 Clips

**C31-10-30-07-E**



Clip Applicator for Polymer Clips  
Dia 10mm, Length 30cm  
(Hem-o-lok Type)

**CH4-10-30-10-PP**

(Purple / Large)


**CH5-10-30-10-PY**

(Yellow / Extra Large)

Clip Applicator for Polymer Clips  
Dia 5mm, Length 30cm  
(Hem-o-lok Type)

**CH3-05-30-10-PG**

(Green / Medium Large)



Clip Applicator for Polymer Clips  
Dia 10mm, Length 30cm (Hem-  
o-lok Type) (Non-Detachable)

**NCH4-10-30-10-PP**

(Purple / Large)

**NCH5-10-30-10-PY**


(Yellow / Extra Large)

Clip Applicator for Polymer Clips  
Dia 5mm, Length 30cm  
(Hem-o-lok Type)

**(Non-Detachable)**

**NCH3-05-30-10-PG**

(Green / Medium Large)




# RETRACTOR




Fan Retractor Take-A-Part  
Dia 5mm, Length 30cm

**FCT-05-30-05**




Fan Retractor Take-A-Part  
Dia 10mm, Length 30cm


**FCT-10-30-10**



# NEEDLE HOLDERS



Axial grip



Pistol Grip



Needle Holder T.C. Straight Jaw  
Axial Handle Dia 5mm, Length 30cm  
■ **NSA-05-30**

Needle Holder T.C. Straight Jaw Pistol  
Handle Dia 5mm, Length 30cm  
■ **NSP-05-30**



Needle Holder T.C. Curve Jaw (Left)  
Axial Handle Dia 5mm, Length 30cm  
■ **NLA - 05 - 30**

Needle Holder T.C. Curve Jaw (Left)  
Pistol Handle Dia 5mm, Length 30cm  
■ **NLP-05-30**



Needle Holder T.C. Tooth Jaw Axial  
Handle Dia 5mm, Length 30cm  
■ **NTA-05-30**

Needle Holder T.C. Tooth Jaw Pistol  
Handle Dia 5mm, Length 30cm  
■ **NTP-05-30**




Needle Holder Self Aligning Jaw  
Axial Handle Dia 5mm, Length 30cm  
■ **NAA-05-30**

Needle Holder Self Aligning Jaw Pistol-  
Handle Dia 5mm, Length 30cm  
■ **NAP-05-30**


**RECOMMENDED Needle SIZE : 30/35mm, 1/2 Circle Round / Taper Cut Needle**

ergonomic designed axial and pistol grip handle gives a good grip in surgeon's hand resulting natural suturing motion


# SUCTION IRRIGATION




2 Way Sliding Valve Suction  
Irrigation with Cannula  
Dia 5mm Length 30 cm  
Shaft Non-Insulated  
■ **SI-05-30-2WC**



2 Way Sliding Valve Suction  
Irrigation with Cannula  
Dia 10mm Length 30 cm  
Shaft Non-Insulated  
■ **SI-10-30-2WC**




Trumpet Valve Suction  
Irrigation with Cannula  
Dia 5mm, Length 30cm  
Shaft Non-Insulated  
■ **SI-05-30-TC**




Trumpet Valve Suction  
Irrigation with Cannula Set  
of 2, Dia 10mm, Length 30cm  
Shaft Non-Insulated  
■ **SI-10-30-TCS**


# KNOT PUSHER




Bhandarkar's knot Pusher  
( $\text{Ø}=1.2\text{mm}$ ) with Reducer  
dia 5mm, Length 40cm  
■ **KPBR-05-40-12**



Bhandarkar's knot Pusher  
( $\text{Ø}=0.9\text{mm}$ ) with Reducer  
dia 5mm, Length 40cm  
■ **KPBR-05-40-09**



Clark's Type Knot  
Pusher ( $\text{Ø}=1.2\text{mm}$ )  
Dia 5mm, Length 40cm  
■ **KPC-05-40-12**



Hook's Type Knot  
Pusher ( $\text{Ø}=1.2\text{mm}$ )  
Dia 5mm, Length 40cm  
■ **KPL-05-40-12**

TYPE OF SUTURE : 1 & 1.0

TYPE OF SUTURE : 2.0



# MONOPOLAR ELECTRODES




'L' Hook Electrode (Teflon)  
Dia 5mm Length 30cm  
■ LHT-05-30




'L' Hook Electrode (Ceramic)  
Dia 5mm Length 30cm  
■ LHC-05-30




Spatula (Teflon)  
Dia 5mm, Length 30cm  
■ SE-05-30



Aspiration Needle Insulated  
Teflon Dia 5mm , Length 30cm  
■ ANIT-05-30-18




Aspiration Needle Non Insulated  
Dia 5mm, Length 30cm  
[Needle Gauge-14]  
■ ANN-05-30-18




Aspiration Needle Non Insulated  
Dia 5mm ,Length 30cm  
[Needle Gauge-16]  
■ ANN-05-30-18


# PORT CLOSURE



Nande's Port Closure  
Left Curve  
■ PCL-12




Nande's Port Closure  
right Curve  
■ PCR-12




Cobbler Needle  
Dia 2mm, Length 11cm  
■ CN-02-11


# SUTURE PASSER



Suture Passer Storz Type  
Dia 2.5mm, Length 11cm  
■ SPT-25-11-M1C




Suture Passer Storz Type  
Dia 2.1mm, Length 11cm  
■ SPT-21-11-M1C



Suture Passer Cone (Teflon)  
SPC-10



# LAP BI CLAMP GRASPER



Lap Bi Clamp Dissector Maryland dia 5mm, length 30cm with Bipolar Plain Handle of PR Series  
 ■ LBMC-05-30-21

Lap Bi-Clamp Plain Grasper (DA) dia 5mm, length 30cm with Bipolar Plain Handle of PR Series  
 ■ LBPC-05-30-21

Lap Bi-Clamp Fenestrated Grasper (DA) dia 5mm, length 30cm with Bipolar Plain Handle of PR Series  
 ■ LBFC-05-30-21

Lap Bi Clamp Curve Partial Front End Sharp Edge Jaw [Cutter] (DA) dia 5mm, Length 30cm with Bipolar Plain Handle of PR Series  
 ■ LCPC-05-30-21

## Gynaec Instruments



Uterine Manipulator Tooth Type Round Colpotomiser with 3 Tips  
 ■ TUM - FC  
 ■ TUM - HC




Vaginal Tube (Small) front end teflon  
 Dia 32mm, Length 20cm  
 ■ VTST-32-20-C



Vaginal Tube (Medium) front end teflon  
 Dia 38mm, Length 20cm  
 ■ VTMT-38-20-C

## Myoma screw




Myoma Screw  
 Dia 5mm, Length 30cm  
 ■ MYS-05-30




Myoma Screw with Needle  
 Dia 5mm, Length 30cm  
 ■ MYN-05-30

## Myomectomy Needle




Myomectomy Needle  
 Dia 5mm, Length 30cm & 17cm  
 ■ MN-05-30-C  
 ■ MN-05-17-C

## Ring applicator



Ring Applicator Storz Type [For DP]  
 Dia 7mm, Length 33cm with Loader Pusher

■ RSD-07-33-C  
 ■ RAL-01 (Loader)  
 ■ RAP-01 (Pusher)



Ring Applicator Storz Type [For SP]  
 Dia 6mm, Length 33cm with Loader Pusher

■ RSS-06-33-C  
 ■ RAL-02 (Loader)  
 ■ RAP-01 (Pusher)

## Endoknife



Endoknife Retractable Detachable Dia 5mm, Length 30cm Complete

■ ER-05-30-C ■ ER-10-30-C  
 ■ EK-05 (Blade) ■ EK-10 (Blade 10/11)

## pcod needle



PCOD Needle New TAP  
 Dia 5mm. Length 30cm  
 ■ PCOD-05-30-D